Manual Runner Handbike (Electro)
Udinumber 870261760Run9700005321SF

Features Runner Handbike

- Both junior and senior.
 - Including adapter (under the wheelchair).
 - For folding wheelchair: adapter detachable.
 - Stainless steel frame.
 - Height and length adjustable.
 - 20" aluminium rim.
 - 7 Gears.
 - 2nd brake/Parking brake.
 - Coaster brake.
 - Ergonomic handles. • Price corrector.
 - Very easy to connect (roller system).
 - Wedge angle and center of gravity of the wheelchair remain adjustable.
 - Bike stand after disconnection
 - Lighting
 - CE marking.
 - Weight: ± 8.5 kg, titanium ± 7.5 kg.
 - Warranty: parts 1 year / frame 5 years.
 - Delivery time: 5 working days after receipt of the wheelchair at Roam.
- Installation directly by appointment is also possible.

ROAM Special Cycles B.V.
Buitenkamp 20
5321KB Hedel
Tel.nr. 073-5997042
E-mail: info@roam.nl
Internet: www.roam.nl

Introduction

Roam is a manufacturer of unique and special "bicycles" whereby independency of users is priority number one.

The Roam Runner Handbike is available in both junior as well as a senior version and can easily be connected to any wheelchair.

Roam was founded in 1987 by Theo de Rouw, who started doing business - on his own initiative - in the evenings. His goal was: independency of disabled people as much as possible.

As a result of this objective the business grew into a full-time job in less than 3 years. With the same objective in mind, the company was converted into a V.O.F. in 1997 and wife Johanna de Rouw became one of the partners. Later on Johan de Rouw joined the company.

Over the years until now independency of disabled people has still remained the most important keyword for Roam.

The management director is Johan de Rouw and the company has 12 staff members. The company is NEN-EN-ISO 9001:2015 and NEN-EN-ISO 13485 certified and works according to the National Care Standards. These certificates ensure that Roam delivers quality - from development and design to delivery and aftercare- .

Roam is a custom friendly company. Her product range is subject to continuous improvement in relation to regularly compliance and user's wishes and needs.

In 2003, Roam was proclaimed: the most creative entrepreneur of the year by the Chamber of Commerce (the certificate was awarded by the Chamber of Commerce). The Roam approach now enjoys international prestige.

Welcome to the Roam-team

Thank you for purchasing our product, the Runner Handbike.

With great care, the Roam team has made your Runner Handbike roadworthy. The user manual has already been received and describes how the Runner hand bike should be handled. Please note that the complete manual should be read. Roam is not responsible for incorrect use and ignoring the commandments / warnings. In case of incorrect or improper use of the hand bike or manual Roam cannot be held responsible. In other words, as a user you must follow the commandment described in the Manual.

The adventure of freedom and a step towards independence can begin with the Runner Handbike and that means: first getting acquainted with this beautiful product; follow the user manual. Take your time to experience!

The Runner Handbike is made to serve you as well as possible and we hope we enjoy your bike. Do not hesitate to make your findings known to your dealer or the manufacturer.

We wish you a lot of fun with the Roam Runner Handbike.

ROAM Special Cycles B.V.
Buitenkamp 20
5321KB Hedel
Tel.nr. 073-5997042
E-mail: info@roam.nl
Internet: www.roam.nl

Index

Introduction 2

Welcome to the Roam Team 3

Table of Contents 4

Control and Delivery 5

Commandments and warnings 6/7/8/9

Sharp adjustment of the wheelchair 9

Connecting and disconnecting the Handbike 9

Brakes 10

Gears 10

Maintenance and your professional 11

Maintenance 11

Tools intended for dealer/Roam Specialist 11

Checklist 12

Tightening moments for dealer/Roam Specialist 12

Warranty 13

EU Declaration of conformity 14

Delivery check

Check that the delivery of the hand bike matches the packing slip. If transport damage is detected upon delivery, please inform your dealer. Roam trusts that the Runner Handbike meets your expectations and wishes you a lot of success and above all a lot of cycling pleasure.

Important

Upon delivery of your Runner Handbike, all parts described on the attached CD are tailor made to your wishes. Certain parts need to be readjusted after some time. The best thing to do is making an appointment with your dealer.

Commandments and warnings

Always observe the commandments and warnings below!

General

1. When purchasing the Runner Handbike, the dealer must provide the correct information to ROAM Special Cycles BV about the use and application.
2. Before you start cycling, always check whether the adapter under the wheelchair is properly secured. In the case of detachable adapters, the star buttons or quick-release fasteners must be checked to see if they are properly secured.
3. Roam cannot be held liable for defects, damage, unsafe actions or situations, accidents or incidents due to incorrect use, incorrect information or incorrect maintenance or repair.
4. Right-angled and sharp turns should be taken at an appropriate speed to avoid the risk of tipping over.

No changes

1. No changes may be made to the Runner. Roam is not responsible for damage caused to the Runner by making changes.
2. The speed of our products may not be increased or changed in any way. When changing our product, the manufacturer is in no way responsible or liable.

User

1. The user of the Runner Handbike should be in good mental health and not be physically restricted in a sitting position.
2. The driver must be able to operate the Runner Handbike properly. The driver's hearing and vision must be suitable for the use of the Runner Handbike. The user should read and understand the manual carefully.
3. If the user is so limited that he/she has vision/hearing problems and/or cannot read and understand the manual, parents/guardians/supervisors are responsible for the application of the manual and the use of the Handbike by the end user.
4. The user is responsible and liable for proper use, as stated in this manual.
5. It is forbidden to use the Runner Handbike under the influence of alcohol / drugs or mind-influencing medication.
6. When biking, make sure you are well equipped (mobile phone).
7. Include your dealer's number in your mobile phone. Provide good further equipment such as: repair-set tires, small first aid kit, and water bottle with drinking water, some ransom and identity card.
8. Provide good further equipment such as: repair kit tires, small first aid kit, and water bottle with drinking water, some ransom and the identity card.
9. When using the Runner Handbike in the dark/dark/twilight or bad weather conditions, lighting should be provided.
10. On longer journeys, protect yourself from bad weather conditions in winter or excessive heat in the summer period.
11. The Runner Handbike is personal. It is forbidden for unauthorized persons to use the Runner Handbike.
12. Only the owner/user is allowed to use the Runner Handbike.

Regulations

1. Adhere to the Bicycle Act which applies in the country.
2. Comply with the Bicycle Act that applies abroad (in the country you live) and adapt to the circumstances.
3. For more information, please also visit websites concerned.
4. Comply with the laws and regulations of the central government.

Biking

1. A user/owner must be at least 16 years of age or older to be able to independently drive the above-mentioned powered products because the speed of these products is more than 10 km per hour.
2. Under the age of 16, the user must be under the supervision of a parent or guardian when the Handbike Runner is used.
3. The user may use: the sidewalk, the footpath, the cycle path/ moped path and the roadway. If all road components are present, the user/owner may decide for himself which road is to be taken.
4. It is forbidden to use the motorway. Follow the traffic rules of the road being used. On the sidewalk, for example, you can drive a maximum of 6 kilometers per hour. Electronic devices, such as mobile phones, navigation systems, tablets or music players, should not be held or used while driving.
5. The same rules apply as when texting, calling or listening to music on the bike. This refers to the website of the Central Government.
6. Biking with the Runner. Adjust your speed to the conditions that occur. Think of busy and bad roads, obstacles and roads along water features.
7. Keep your hands on the wheel and don't get distracted.
8. Before shifting gears while driving, the pedals must be kept still.
9. Take obstacles in the road surface (thresholds) at an appropriate speed.
10. The hands should always be held on the handles during biking (with the exception of shifting gears).
11. Do not put the hands in rotating parts.
12. For good visibility in traffic, a flag is flown.
13. To protect the head, a helmet is worn.
14. For maintenance, settings, adjustments, replacement of parts and / or repairs, reference is made to the dealer / Roam specialist.
15. The Roam specialist as well as the dealer is obliged to follow the manual and to apply the recommended tools.
16. For roads where extra attention is required such as roads: with obstacles, along waterways, equipped with elevations, railway crossings, tram rails or pedestrian crossing and or driving along busy roads etc., must be driven at walking pace.
17. Driving should always be done at an appropriate speed, depending on the circumstances.

Slopes and brakes

1. In a number of cases, if this manual is not followed correctly, a dangerous situation may arise for the user and passengers of the Runner with regard to taking slopes.
2. For slopes of 5% or more, the brakes are not calculated. On slopes that are very long or steep, the brakes can run hot and the engine can become overloaded so that defects occur.
3. The preceding point regarding slopes also applies to descents. Note that when descending a slope, the speed increases.
4. Note: On roads where oil, mud, sand or leaves are on the road surface, the braking distance may be longer. This also applies to uneven road surfaces (potholes/holes) and winter conditions.

Curves/Angles

Right-angled and sharp turns should be taken at an appropriate speed (depending on the circumstances: road and physical capabilities) to avoid the risk of tipping over.

Battery

1. Store the battery at a temperature below 25°C. The battery must be at least 30% and a maximum of 80% charged.
2. The battery must be charged in the fireproof case, which is purchased from the Runner.
3. The battery, when not in use, is stored and stored in the fireproof case.
4. Make sure that the battery is not discharged below 30% when in use.
5. Use only original chargers (no fast chargers).
6. Keep an eye on the battery while charging; never charge at night due to fire risk in case of overheating. If a cell of the lithium battery catches fire, the other cells also catch fire.
7. A lithium battery is extremely difficult to extinguish. In case of fire, call the alarm number immediately.
8. A faulty battery can leak and release fumes. Sparks can also occur or the battery can smell strong and become extremely hot. In one of these cases, disconnect the battery.
9. If a device with a lithium ion battery leaks, gets excessively hot, smells strong, produces smoke or produces sparks, it is wrong. If you notice any of the signals mentioned, it is important that you disconnect the device from the charger and move it outside.
10. Make sure you do not place the device near flammable materials and then call the alarm number.
11. Replace damaged chargers.
12. Do not store batteries in direct sunlight.
13. When purchasing the Runner, the dealer must provide the correct information to ROAM Special Cycles BV about the use and application.

Environment

A defective, leaking, fallen, or deformed battery battery should be replaced immediately. The old battery must be returned to the dealer or to the environmental street of the municipality.

Sharp adjustment of the wheelchair

Warning: Pay attention when adjusting the wheelchair sharply. If your wheelchair is sharply adjusted, the Runner Handbike does not have as much wheel pressure and the braking distance is many times longer. The user must take this into account!

Handbike

Connecting

Connecting the Runner Handbike under your wheelchair is very easy. First of all, you put the wheelchair on the brakes; then you close the lever of the connector on the Runner Handbike (down) and now slide the Runner Handbike with the bearing holder as far as possible into the adapter of your wheelchair. Now you need to open the lever of the connector on the Runner Handbike (up), so that the bracket falls over the tab of the adapter, then close the trade again. Now the Runner Handbike is docked. Put the wheelchair off the brake and you can start driving. If there are any ambiguities for you, we refer you to the videos on our website: www.:Roam.nl

Disconnecting

Disconnecting the Runner Handbike is also very easy. First put the wheelchair on the brake, and then you have to stop the Runner Handbike, preferably grip the handlebars, and then open the connector with your other hand and hold up the bracket of the clamp. Now you can close the lever of the connector again so that the bracket is on the tab. Then slowly slide the Runner Handbike out of the adapter and put it away. If there are any ambiguities for you, we refer you to the videos on our website: www.:Roam.nl

Brakes

It is vital that you fully understand the operation of your Runner Hand bike's braking system. If the braking system is not used correctly, this may result in you losing control of the steering wheel or may lead to serious consequences (accident, serious injuries). The Runner Hand bike is equipped with a disc brake. If, in the opinion of the user, the brake does not function properly, the dealer should be contacted immediately. The dealer must resolve any defects. It is forbidden to use the Runner Handbike when the braking system does not work in any form.

Note: The disc brake can get hot; contact the dealer immediately if there is any doubt as to whether the brake is not working correctly. Do not come into contact with the disc brake, construction or housing. Before biking, always check your brake.

Coaster brakes

The drive is done by turning the pedals forward. The brake is activated by turning the pedals in the opposite direction: so turn back! Pay attention to slopes: if these are more than 5%, the brake can run hot and no longer brake! In this situation, you must keep the speed limited (no more than 20 km per hour).

Gears

The bike comes with 7 Shimano gear. The gears can only be used while driving. For shifting while driving, one should keep the pedals still Shift with one gear at a time and reduce during the shifting force you apply to the pedals. If you force-turn the pedals around, trying to force the shift lever to operate, there is a chance that you will slide off the pedals, the bike will overturn and you will be seriously injured.

Maintenance and your professional

Vital parts such as cables, spokes and brakes get a chance to "set" themselves in the start-up period. However, in order to be able to continue driving pleasantly and safely, a major annual turn is necessary. Your professional not only has the tools and experience to set and lubricate hidden rotating parts such as headset and hubs, but his hands, eyes and ears are also able to detect small imperfections in time. Before winter sets in or before a holiday trek, we advise you to have the bike checked. Proper maintenance is also extremely important with regard to the warranty. Make an appointment with your dealer before the service will be executed.

Maintenance Runner Handbike

1. Clean the Runner Handbike with a dry cloth.
2. In case of serious dirt, clean the Runner Hand bike with a soft damp sponge and dry it.
3. Do not spray the Runner Handbike with a high pressure washer.
4. Do not apply chemical cleaning agents (chlorine/bleach/soda) when cleaning the Runner Handbike.
5. Lubricate aluminium and chrome parts with an acid-free petroleum jelly or special grease (ask your dealer).
6. Store the Runner Handbike dry.

Tools intended for dealer/Roam Specialist

For the modest work, a simple set of tools is sufficient. You are already well on your way with a set of wrenches or socket wrenches (8 to 15), Allen wrenches, screwdriver, combination pliers, a can of bicycle oil and of course a tire repair kit. When replacing, make sure that original or equivalent parts are assembled for better operation and because in some cases the warranty will otherwise expire. Ask your Roam dealer for advice!

Checklist

Check regularly:

1. Are the brakes working properly?
2. Are the brake pads not worn out?
3. Is stem not higher than the safety marking indicates?
4. Are all the nuts and bolts safe?
5. Is the tread of the tires not worn out and do they not have weak spots or bumps?
Are the tires on pressure (tire pressure is on the side of the tire)?
6. Is the stem and steering rod tightly attached?
7. Is the lighting working properly and are the batteries of the rear light not due for replacement?
8. Is the side reflection of the tires clean and therefore clearly visible in the dark? Are the handles still secure?
9. Are the reflectors of the pedals free of dirt and are the pedals still sufficiently stiff?
10. Is the quick release working properly? (there must be no play on it, nor must the quick release be loose).
11. If in doubt, contact your dealer.

Tightening moments intended for dealer/Roam Specialist

Below you will find an overview of recommended tightening moments. We recommend that you tighten the listed parts with a torque wrench based on the values mentioned.

Tightening the parts harder can lead to damage to your bike. It is best to ask your Roam dealer for help and advice for this action.

Steering plug bolt 12-15 Nm

Expander bolt with angled bobbin 20-23 Nm

Hub axle nuts front wheels 20-27 Nm

Guarantee

Roam gives a 5-year warranty against defects in material and finish on the frame. All parts and components manufactured and marketed by Roam have a 1-year warranty. This warranty does not include: tire and parts wear. Roam will repair or replace parts where necessary, provided that an official Roam dealer has agreed to return within the established warranty period. Roam does not guarantee parts that are damaged as a result of incorrect assembly and / or repair by the relevant owner, neglect, accident or incorrect use.

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5321KB Heden

E-mail: info@roam.nl.

SRN (Single Registration Number): APP000006468

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Declares that the below stated Roam Product:

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English

1. **an active device that provides relief, improvement or support for physical immobility, disabilities or limitations.**
2. **expanding/improving independent functioning as much as possible.**

Udi-nummer (Unique identification number): **870261760Run97000005321SF**

Risicoklasse: **1** (conform de classificatie MDR 2017/745 EU).

Risk classification **1**: (according classification MDR 2017/745 EU).

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Comply with the below stated European regulation.

- **MDR (Medical Device Regulation)(2017/745/ EU)**
- **Machinerichtlijn(2006/42/EG) - Machinery Directive**
- **EMC richtlijn(2014/30/EU) - Electric Magnetic Compatibility Directive**
- **Laagspanningsrichtlijn(2014/35/EU)- Low voltage Directive**

Datum (Date): **1 January 2023**

Handtekening (Signature):

Roam Special Cycles

Johan de Rouw, Managing Director/JDR